

FITC-labeled ODN 1585 sodium

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	FITC-labeled ODN 1585 (sodium) effectively induces IFN and TNF α production and serves to assess CpG ODN cellular uptake and localization via confocal laser-scanning microscopy [excitation 495 nm, emission 520 nm] or flow cytometry.
Targets(IC50)	Others,IFNAR,TLR

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481