

Biotin-labeled ODN 1668 sodium

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Biotin-labeled ODN 1668 (sodium), a class B CpG oligodeoxynucleotide (ODN) and TLR9 agonist, enables the evaluation of cellular uptake and localization of CpG ODNs using a biotin detection system and light microscopy.
Targets(IC50)	Others,TLR

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481