

## Biotin-labeled ODN 2216 sodium

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Biotin-labeled ODN 2216 (sodium) acts as a specific agonist for human TLR9 (toll-like receptor 9), enabling the evaluation of CpG ODN cellular uptake and localization via a biotin detection system coupled with light microscopy.
Targets(IC50)	Others,TLR

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481