

FITC-labeled ODN 2216 sodium

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	FITC-labeled ODN 2216 (sodium), a specific agonist for human toll-like receptor 9 (TLR9), facilitates the assessment of CpG ODN cellular uptake and localization through confocal laser-scanning microscopy (excitation 495 nm, emission 520 nm) or flow cytometry.
Targets(IC50)	Others,Interleukin,TLR

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481