

FITC-labeled ODN TTAGGG sodium

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	FITC-labeled ODN TTAGGG (sodium) is an inhibitory oligonucleotide that antagonizes TLR9, AIM2, and cGAS, aiding in assessing CpG ODN cellular uptake and localization through confocal laser-scanning microscopy (excitation at 495 nm, emission at 520 nm) or flow cytometry.
Targets(IC50)	Others,TLR

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481