


## VEGFR2-IN-3

## Chemical Properties

CAS No. : 417717-09-0  
 Formula: C<sub>26</sub>H<sub>28</sub>ClN<sub>5</sub>O<sub>4</sub>  
 Molecular Weight: 509.98  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	VEGFR2-IN-3 (compound 385) is a potent VEGFR2 inhibitor [1].
Targets(IC50)	VEGFR

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9609 mL	9.8043 mL	19.6086 mL
5 mM	0.3922 mL	1.9609 mL	3.9217 mL
10 mM	0.1961 mL	0.9804 mL	1.9609 mL
50 mM	0.0392 mL	0.1961 mL	0.3922 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481