


TRPC3/6-IN-2

## Chemical Properties

CAS No. : 2387893-55-0  
 Formula: C<sub>18</sub>H<sub>23</sub>F<sub>2</sub>N<sub>5</sub>  
 Molecular Weight: 347.41  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	TRPC3/6-IN-2 is a potent inhibitor of TRPC3 and TRPC6, with IC <sub>50</sub> values of 16 nM and 29.8 nM, respectively [1].
Targets(IC <sub>50</sub> )	Others,TRP/TRPV Channel

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8784 mL	14.3922 mL	28.7844 mL
5 mM	0.5757 mL	2.8784 mL	5.7569 mL
10 mM	0.2878 mL	1.4392 mL	2.8784 mL
50 mM	0.0576 mL	0.2878 mL	0.5757 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481