

## PROTAC SOS1 degrader-4

## Chemical Properties

CAS No. :

Formula: C<sub>53</sub>H<sub>65</sub>ClN<sub>8</sub>O<sub>7</sub>S

Molecular Weight: 993.65

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	PROTAC SOS1 Degrader-4 (compound 10) is a potent degrader of SOS1 and exhibits antiproliferative activity [1].
Targets(IC50)	PROTACs,Ras

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0064 mL	5.032 mL	10.0639 mL
5 mM	0.2013 mL	1.0064 mL	2.0128 mL
10 mM	0.1006 mL	0.5032 mL	1.0064 mL
50 mM	0.0201 mL	0.1006 mL	0.2013 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481