

CDK9-IN-24

Chemical Properties

CAS No. :

Formula: C₂₇H₃₁N₃O₄

Molecular Weight: 461.55

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	CDK9-IN-24 (compound 21a) is a potent and selective CDK9 inhibitor that significantly inhibits tumor proliferation, impedes cell growth, and triggers apoptosis via downregulation of Mcl-1 and c-Myc, making it valuable for acute myeloid leukemia research [1].
Targets(IC50)	Apoptosis,CDK

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1666 mL	10.8331 mL	21.6661 mL
5 mM	0.4333 mL	2.1666 mL	4.3332 mL
10 mM	0.2167 mL	1.0833 mL	2.1666 mL
50 mM	0.0433 mL	0.2167 mL	0.4333 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481