

FAK-IN-11

Chemical Properties

CAS No. :

Formula: C₂₅H₄₁N₃O₃

Molecular Weight: 403.6

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	FAK-IN-11 (Compound 4l), a FAK inhibitor, specifically targets the ATP binding pocket of FAK to prevent its phosphorylation. It exhibits cytotoxic effects on MDA-MB-231 cells, achieving an IC ₅₀ value of 13.73 μ M, and induces non-apoptotic cell death in these cells [1].
Targets(IC ₅₀)	FAK

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4777 mL	12.3885 mL	24.777 mL
5 mM	0.4955 mL	2.4777 mL	4.9554 mL
10 mM	0.2478 mL	1.2389 mL	2.4777 mL
50 mM	0.0496 mL	0.2478 mL	0.4955 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481