


## Nevadistinel

## Chemical Properties

CAS No. : 2181816-92-0  
 Formula: C<sub>15</sub>H<sub>18</sub>N<sub>2</sub>O<sub>3</sub>  
 Molecular Weight: 274.32  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Nevadistinel (NYX-458; NYX-3054) is a positive allosteric modulator of the N-methyl-D-aspartate (NMDA) receptor, indicated for mitigating cognitive deficits related to neurodegenerative conditions, including mild cognitive impairment, mild Alzheimer's disease, Parkinson's disease, and Lewy body disease [1].
Targets(IC50)	iGluR

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6454 mL	18.2269 mL	36.4538 mL
5 mM	0.7291 mL	3.6454 mL	7.2908 mL
10 mM	0.3645 mL	1.8227 mL	3.6454 mL
50 mM	0.0729 mL	0.3645 mL	0.7291 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481