


## Risevistinel

## Chemical Properties

CAS No. :	2591344-26-0
Formula:	C <sub>14</sub> H <sub>23</sub> N <sub>3</sub> O <sub>4</sub>
Molecular Weight:	297.35
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Risevistinel (NYX-783), a positive allosteric modulator of the N-methyl-D-aspartate (NMDA) receptor, may mitigate cognitive deficits in neurodegenerative conditions, including mild cognitive impairment, mild Alzheimer's disease, Parkinson's disease, and Lewy body disease [1].
Targets(IC50)	iGluR

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.363 mL	16.8152 mL	33.6304 mL
5 mM	0.6726 mL	3.363 mL	6.7261 mL
10 mM	0.3363 mL	1.6815 mL	3.363 mL
50 mM	0.0673 mL	0.3363 mL	0.6726 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481