


GDC-6599

## Chemical Properties

CAS No. : 2376824-99-4  
 Formula: C<sub>20</sub>H<sub>19</sub>ClN<sub>6</sub>O<sub>3</sub>  
 Molecular Weight: 426.86  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
*Actual storage temperature shall be subject to the COA.*



## Biological Description

Description	GDC-6599 (Example 8) is an orally active TRPA1 channel inhibitor, potentially useful for researching TRPA1-mediated conditions, including pain [1].
Targets(IC50)	Others,TRP/TRPV Channel

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3427 mL	11.7134 mL	23.4269 mL
5 mM	0.4685 mL	2.3427 mL	4.6854 mL
10 mM	0.2343 mL	1.1713 mL	2.3427 mL
50 mM	0.0469 mL	0.2343 mL	0.4685 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481