

Bak BH3 (72-87), TAMRA-labeled

Chemical Properties

CAS No. :

Formula: C97H145N27O28

Molecular Weight: 2137.35

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Bak BH3 (72-87), TAMRA-labeled, is a bioactive peptide.
Targets(IC50)	Bcl-2 Family

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.4679 mL	2.3393 mL	4.6787 mL
5 mM	0.0936 mL	0.4679 mL	0.9357 mL
10 mM	0.0468 mL	0.2339 mL	0.4679 mL
50 mM	0.0094 mL	0.0468 mL	0.0936 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481