


Bim BH3

Chemical Properties

CAS No. :	505070-06-4
Formula:	C108H170N32O31S
Molecular Weight:	2444.77
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Bim BH3 is a bioactive peptide that belongs to the pro-apoptotic Bcl-2 family of proteins.
Targets(IC50)	Bcl-2 Family

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.409 mL	2.0452 mL	4.0904 mL
5 mM	0.0818 mL	0.409 mL	0.8181 mL
10 mM	0.0409 mL	0.2045 mL	0.409 mL
50 mM	0.0082 mL	0.0409 mL	0.0818 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481