

Beta-glucuronidase (bovine liver)

Chemical Properties

CAS No. : 9001-45-0

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Beta-glucuronidase, a critical lysosomal enzyme, is implicated in the breakdown of glucuronate-containing glycosaminoglycans [1].
Targets(IC50)	Others
In vitro	Beta-glucuronidase deficiency leads to mucopolysaccharidosis type VII (MPSVII), resulting in lysosomal accumulation within the brain [1].

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481