


Tyr-(D-Dab4,Arg5,D-Trp8)-cyclo-Somatostatin-14 (4-11)

Chemical Properties

CAS No. :	496849-46-8
Formula:	C67H85N15O11
Molecular Weight:	1276.49
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Tyr-(D-Dab4,Arg5,D-Trp8)-cyclo-Somatostatin-14 (4-11) is a somatostatin agonist used in cancer research [1].
Targets(IC50)	Somatostatin

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.7834 mL	3.917 mL	7.834 mL
5 mM	0.1567 mL	0.7834 mL	1.5668 mL
10 mM	0.0783 mL	0.3917 mL	0.7834 mL
50 mM	0.0157 mL	0.0783 mL	0.1567 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481