

## Tripeptide-41

## Chemical Properties

CAS No. : 1093241-16-7


Formula: C<sub>29</sub>H<sub>30</sub>N<sub>4</sub>O<sub>5</sub>

Molecular Weight: 514.57

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Tripeptide-41 (CG-Lipoxyn) is a bioactive peptide known for its ability to reduce fat accumulation. It has been utilized as an ingredient in cosmetic formulations [1].
Targets(IC50)	NF-κB

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9434 mL	9.7169 mL	19.4337 mL
5 mM	0.3887 mL	1.9434 mL	3.8867 mL
10 mM	0.1943 mL	0.9717 mL	1.9434 mL
50 mM	0.0389 mL	0.1943 mL	0.3887 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481