

Bid BH3 (80-99), FAM labeled

## Chemical Properties

CAS No. :

Formula: C116H171N33O38S

Molecular Weight: 2667.86

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	FAM-labeled Bid BH3 (80-99) is a biologically active peptide featuring a 5-FAM label at its N-terminus.
Targets(IC50)	Bcl-2 Family

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.3748 mL	1.8742 mL	3.7483 mL
5 mM	0.075 mL	0.3748 mL	0.7497 mL
10 mM	0.0375 mL	0.1874 mL	0.3748 mL
50 mM	0.0075 mL	0.0375 mL	0.075 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481