


## Noxa A BH3

## Chemical Properties

CAS No. :	505070-09-7
Formula:	C102H162N26O29S
Molecular Weight:	2248.6
Storage:	Keep away from moisture Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	Noxa A BH3, a bioactive peptide, encompasses the BH3 domain of the Noxa protein.
Targets(IC50)	NADPH

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.4447 mL	2.2236 mL	4.4472 mL
5 mM	0.0889 mL	0.4447 mL	0.8894 mL
10 mM	0.0445 mL	0.2224 mL	0.4447 mL
50 mM	0.0089 mL	0.0445 mL	0.0889 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481