


Noxa B BH3

Chemical Properties

CAS No. :	505070-10-0
Formula:	C95H164N30O31S
Molecular Weight:	2254.57
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Noxa B BH3, a biologically active peptide encompassing the BH3 domain of the Noxa protein, displays targeted biological activity.
Targets(IC50)	NADPH

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.4435 mL	2.2177 mL	4.4354 mL
5 mM	0.0887 mL	0.4435 mL	0.8871 mL
10 mM	0.0444 mL	0.2218 mL	0.4435 mL
50 mM	0.0089 mL	0.0444 mL	0.0887 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481