

Hrk BH3

Chemical Properties

CAS No. : 1644272-97-8


Formula: C99H160N30O31

Molecular Weight: 2266.51

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Hrk BH3, a biologically active peptide comprising the BH3 domain of Hrk, exhibits functionality within cellular processes.
Targets(IC50)	Bcl-2 Family

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.4412 mL	2.206 mL	4.4121 mL
5 mM	0.0882 mL	0.4412 mL	0.8824 mL
10 mM	0.0441 mL	0.2206 mL	0.4412 mL
50 mM	0.0088 mL	0.0441 mL	0.0882 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481