


Aurein 2.4

Chemical Properties

CAS No. :	302343-07-3
Formula:	C77H133N19O19
Molecular Weight:	1629
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Aurein 2.4, an antibiotic antimicrobial peptide, exhibits activity against Gram-positive bacteria [1] [2].
Targets(IC50)	Antibacterial, Antibiotic

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.6139 mL	3.0694 mL	6.1387 mL
5 mM	0.1228 mL	0.6139 mL	1.2277 mL
10 mM	0.0614 mL	0.3069 mL	0.6139 mL
50 mM	0.0123 mL	0.0614 mL	0.1228 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481