


## Bombinin H5

## Chemical Properties

CAS No. :	325686-32-6
Formula:	C91H165N23O21
Molecular Weight:	1917.43
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Bombinin H5, an antimicrobial peptide originating from the skin of the (moth) <i>Bombina variegata</i> [1], exhibits potent biocidal properties.
Targets(IC50)	Others,Antibacterial

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.5215 mL	2.6077 mL	5.2153 mL
5 mM	0.1043 mL	0.5215 mL	1.0431 mL
10 mM	0.0522 mL	0.2608 mL	0.5215 mL
50 mM	0.0104 mL	0.0522 mL	0.1043 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481