

## Anti-Mouse CD117 Antibody

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from moisture  
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Anti-Mouse CD117 Antibody, a rat-derived IgG2b isotype, acts as an inhibitor specific to the mouse CD117 antigen.
Targets(IC50)	c-Kit

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481