

Eurestobart

Chemical Properties

CAS No. : 2682847-53-4

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Eurestobart is a humanized IgG1k (IgG1 kappa) monoclonal antibody that specifically targets the enzyme ectonucleoside triphosphate diphosphohydrolase (ENTPDase).
Targets(IC50)	Phosphatase

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481