

Anti-Mouse 4-1BB Antibody (3H3)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from moisture
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-Mouse 4-1BB Antibody is a rat-derived IgG2a isotype inhibitor targeting the mouse 4-1BB (CD137) protein.
Targets(IC50)	TNF

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481