

Anti-Mouse Fas Ligand Antibody (MFL3)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from moisture

-20°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-Mouse Fas Ligand Antibody, an Armenian hamster-derived IgG inhibitor, targets the mouse Fas Ligand.
Targets(IC50)	Fatty Acid Synthase

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481