

Ledelabricin alfa

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Ledelabricin alfa, an isoform A variant of human proteoglycan 4 (PRG4), functions as a lubricant for joints and boundaries. It is utilized in the study of arthropathy [1].
Targets(IC50)	Endogenous Metabolite

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481