

## Peginterferon lambda-1a

### Chemical Properties

CAS No. : 914617-98-4

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Peginterferon lambda-1a (BMS-914143) is a monoclonal antibody targeting IFNLR1, composed of IFNL1 conjugated to pegol [1].
Targets(IC50)	IFNAR

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481