

Pegsunercept

Chemical Properties

CAS No. : 330988-75-5

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Pegsunercept (PEG sTNF-RI), a pegylated monoclonal antibody, selectively binds to TNFA, incorporating a polyethylene glycol (pegol) moiety [1].
Targets(IC50)	TNF

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481