

Mecbotamab

Chemical Properties

CAS No. : 2460400-28-4

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Mecbotamab vedotin (BA3011) is a humanized IgG1-κ antibody engineered to target the AXL receptor tyrosine kinase (AXL). Serving as a conditionally active biologic (CAB), mecbotamab is conjugated to monomethyl auristatin E (MMAE) through a cleavable linker, creating an antibody-drug conjugate (ADC) for therapeutic use [1] [2].
Targets(IC50)	ADC Antibody,TAM Receptor

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481