

Romilkimab

Chemical Properties

CAS No. : 1399584-78-1

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Romilkimab (SAR156597) is a chimeric humanized immunoglobulin (Ig) antibody specifically designed to target interleukins 4 and 13 (IL-4 and IL-13) [1].
Targets(IC50)	Interleukin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481