

Sasanlimab

Chemical Properties

CAS No. : 2206792-50-7

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Sasanlimab represents a highly specialized humanized IgG4 monoclonal antibody specifically engineered to target the human programmed cell death protein 1 (PD-1) with exceptional binding affinity. It effectively facilitates the measurable blockade of the PD-1/PD-L1 and PD-1/PD-L2 inhibitory axes across various preclinical experimental models to evaluate the restoration of T-cell mediated anti-tumor immunity and the reversal of immune exhaustion during strictly monitored laboratory observation periods.
Targets(IC50)	PD-1/PD-L1

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481