


Androgen receptor degrader-1

Chemical Properties

CAS No. :	2616553-35-4
Formula:	C ₁₅ H ₁₄ ClN ₃ O ₄
Molecular Weight:	335.74
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Androgen Receptor Degrader-1 (Compound 18) is a potent agent for androgen receptor degradation, applicable in cancer research [1].
Targets(IC50)	Androgen Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9785 mL	14.8925 mL	29.785 mL
5 mM	0.5957 mL	2.9785 mL	5.957 mL
10 mM	0.2978 mL	1.4892 mL	2.9785 mL
50 mM	0.0596 mL	0.2978 mL	0.5957 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481