


Dehydroborapetoside B

Chemical Properties

CAS No. :	1221178-16-0
Formula:	C ₂₇ H ₃₄ O ₁₂
Molecular Weight:	550.55
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Dehydroborapetoside B is a natural product for research related to life sciences. The catalog number is TN3795 and the CAS number is 1221178-16-0.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8164 mL	9.0818 mL	18.1637 mL
5 mM	0.3633 mL	1.8164 mL	3.6327 mL
10 mM	0.1816 mL	0.9082 mL	1.8164 mL
50 mM	0.0363 mL	0.1816 mL	0.3633 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481