


(Rac)-Apremilast-D5

Chemical Properties

CAS No. :	1258597-61-3
Formula:	C ₂₂ H ₂₄ N ₂ O ₇ S
Molecular Weight:	465.53
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	(Rac)-Apremilast-D5 is a deuterium-labeled version of the enantiomer (R)-Apremilast, also known as (R)-CC-10004, which itself is one specific form of Apremilast.
Targets(IC50)	Apoptosis,PDE,TNF

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1481 mL	10.7404 mL	21.4809 mL
5 mM	0.4296 mL	2.1481 mL	4.2962 mL
10 mM	0.2148 mL	1.074 mL	2.1481 mL
50 mM	0.043 mL	0.2148 mL	0.4296 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Synthesis of Substituted β -Functionalised Styrenes by Microwave-Assisted Olefin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481