

$\alpha$ -Amylase-IN-4

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	$\alpha$ -Amylase-IN-4 (Compd 10y) exhibits the highest amylase inhibition with an IC50 value of $17.83 \pm 0.14 \mu\text{g/mL}$ [1].
Targets(IC50)	Others,Amylase

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481