


Z-Pro-Ala

Chemical Properties

CAS No. : 14030-00-3
 Formula: C₁₆H₂₀N₂O₅
 Molecular Weight: 320.34
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Z-Pro-Ala, an acid carboxypeptidase, can be isolated from the grains and leaves of <i>Triticum aestivum</i> L. (wheat) [1].
Targets(IC50)	Others,Lipid

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1217 mL	15.6084 mL	31.2168 mL
5 mM	0.6243 mL	3.1217 mL	6.2434 mL
10 mM	0.3122 mL	1.5608 mL	3.1217 mL
50 mM	0.0624 mL	0.3122 mL	0.6243 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481