

ZH8667

Chemical Properties


CAS No. : 910387-04-1

Formula: C₁₄H₁₄FN

Molecular Weight: 215.27

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	ZH8667 is a TAAR1-Gs agonist that holds potential for schizophrenia research [1].
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.6453 mL	23.2266 mL	46.4533 mL
5 mM	0.9291 mL	4.6453 mL	9.2907 mL
10 mM	0.4645 mL	2.3227 mL	4.6453 mL
50 mM	0.0929 mL	0.4645 mL	0.9291 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481