


## Z-Gly-Gly-Arg-AFC

## Chemical Properties

CAS No. :	76409-80-8
Formula:	C <sub>28</sub> H <sub>30</sub> F <sub>3</sub> N <sub>7</sub> O <sub>7</sub>
Molecular Weight:	633.58
Storage:	Keep away from moisture Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	Z-Gly-Gly-Arg-AFC, also known as a plasminogen activator acrosine substrate with fluorescence, is a biologically active peptide.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5783 mL	7.8917 mL	15.7833 mL
5 mM	0.3157 mL	1.5783 mL	3.1567 mL
10 mM	0.1578 mL	0.7892 mL	1.5783 mL
50 mM	0.0316 mL	0.1578 mL	0.3157 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481