

YYRADD

Chemical Properties

CAS No. : 1082226-36-5


Formula: C₃₄H₄₇N₉O₁₁

Molecular Weight: 757.79

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	YYRADD, a polypeptide, enhances collagen production within the skin's extracellular matrix, thus improving skin appearance and aiding the repair of damaged skin [1].
Targets(IC50)	MMP

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3196 mL	6.5981 mL	13.1963 mL
5 mM	0.2639 mL	1.3196 mL	2.6393 mL
10 mM	0.132 mL	0.6598 mL	1.3196 mL
50 mM	0.0264 mL	0.132 mL	0.2639 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481