


YQEAFRRFFGPV

Chemical Properties

CAS No. : 129356-77-0
 Formula: C73H101N19O17
 Molecular Weight: 1516.7
 Storage: Keep away from moisture
 Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	YQEAFRRFFGPV is a dodecapeptide known for its emulsification properties [1].
-------------	------------------------------------------------------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.6593 mL	3.2966 mL	6.5933 mL
5 mM	0.1319 mL	0.6593 mL	1.3187 mL
10 mM	0.0659 mL	0.3297 mL	0.6593 mL
50 mM	0.0132 mL	0.0659 mL	0.1319 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481