


Y4R agonist-2

Chemical Properties

CAS No. :	2757328-63-3
Formula:	C53H81N19O10
Molecular Weight:	1144.33
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Y4R Agonist-2 is a potent agonist of the Y4 receptor, exhibiting a K_i value of 0.033 nM.
Targets(IC50)	Neuropeptide Y Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.8739 mL	4.3694 mL	8.7387 mL
5 mM	0.1748 mL	0.8739 mL	1.7477 mL
10 mM	0.0874 mL	0.4369 mL	0.8739 mL
50 mM	0.0175 mL	0.0874 mL	0.1748 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481