

XQ2B

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	XQ2B is a cGAS inhibitor that impedes protein-DNA interactions and phase separation. It suppresses antiviral immune responses induced by herpes simplex virus 1 (HSV-1) and augments HSV-1 infection in both in vitro and in vivo settings [1].
Targets(IC50)	Autophagy,cGAS

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481