

WAY-659590

Chemical Properties

CAS No. : 867283-39-4


Formula: C₁₇H₁₆N₂O₄

Molecular Weight: 312.32

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	WAY-659590 is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: 80 mg/mL (256.15 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2018 mL	16.0092 mL	32.0184 mL
5 mM	0.6404 mL	3.2018 mL	6.4037 mL
10 mM	0.3202 mL	1.6009 mL	3.2018 mL
50 mM	0.064 mL	0.3202 mL	0.6404 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Joseph M. Sheridan, et al. CASEIN KINASE 1δ (CK1δ) INHIBITORS. WO2012080729A2.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481