

WAY-648936

Chemical Properties

CAS No. : 796888-73-8


Formula: C₁₇H₂₁N₃O₄

Molecular Weight: 331.37

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	WAY-648936 is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: 13 mg/mL (39.23 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0178 mL	15.0889 mL	30.1777 mL
5 mM	0.6036 mL	3.0178 mL	6.0355 mL
10 mM	0.3018 mL	1.5089 mL	3.0178 mL
50 mM	0.0604 mL	0.3018 mL	0.6036 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Bryson HM, et al. Cefetamet pivoxil. A review of its antibacterial activity, pharmacokinetic properties and therapeutic use. *Drugs*. 1993;45(4):589-621.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481