


WAY-621924

## Chemical Properties

CAS No. : 745789-70-2  
Formula: C<sub>17</sub>H<sub>14</sub>N<sub>2</sub>O<sub>2</sub>S  
Molecular Weight: 310.37  
Storage: Keep away from direct sunlight  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.



## Biological Description

Description	WAY-621924 is a bioactive experimental molecule primarily used for biochemical reaction research. It is frequently used as a pharmacological tool in in vitro experimental models.
Targets(IC50)	Others

## Solubility Information

Solubility	DMSO: 80 mg/mL (257.76 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.222 mL	16.1098 mL	32.2196 mL
5 mM	0.6444 mL	3.222 mL	6.4439 mL
10 mM	0.3222 mL	1.611 mL	3.222 mL
50 mM	0.0644 mL	0.3222 mL	0.6444 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

Jiang H, Lan N, Ma W, et al. Synthesis and evaluation of the antitumor activity of 2-amino-4-tetrahydroindazole-substituted benzamide derivatives as HSP90 inhibitors[J]. Journal of Molecular Structure, 2024, 1300: 137266.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481