

WAY-616296

Chemical Properties

CAS No. : 503065-67-6


Formula: C₁₇H₁₄N₂O₂S

Molecular Weight: 310.37

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	WAY-616296 has potential anticholesterol activity and can be used in the study of diabetes, hypertension, psoriasis and obesity.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.222 mL	16.1098 mL	32.2196 mL
5 mM	0.6444 mL	3.222 mL	6.4439 mL
10 mM	0.3222 mL	1.611 mL	3.222 mL
50 mM	0.0644 mL	0.3222 mL	0.6444 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Hae Young Chung, et al. Novel compound having skin-whitening, anti-oxidizing and ppar activities and medical use thereof: US20160102065A1.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481