

WAY-608119

## Chemical Properties


CAS No. : 50362-93-1

Formula: C<sub>13</sub>H<sub>9</sub>N<sub>3</sub>O

Molecular Weight: 223.23

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	WAY-608119 is a bioactive experimental molecule primarily used for biochemical reaction research. It is frequently used as a pharmacological tool in in vitro experimental models.
Targets(IC50)	Others

## Solubility Information

Solubility	DMSO: 16 mg/mL (71.67 mM), Sonication is recommended. ( < 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.4797 mL	22.3984 mL	44.7968 mL
5 mM	0.8959 mL	4.4797 mL	8.9594 mL
10 mM	0.448 mL	2.2398 mL	4.4797 mL
50 mM	0.0896 mL	0.448 mL	0.8959 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

Mashhoon N, et al. Selective inhibitors of bacterial DNA adenine methyltransferases. J Biomol Screen. 2006 Aug;11(5):497-510.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481